

SynAT.Vis - Quick Start

Check a designed gene for problems in Chlamydomonas. No coding needed.

1

Unzip the folder

Right-click the SynAT.Vis zip and choose "Extract All". Open the SynAT.Vis folder.



2

Double-click RUN_ME to check it works


Windows: RUN_ME.bat Mac: RUN_ME.command
A window opens. Look for the green line:

```
> running self-test
6/6 checks passed
the tool works.
```

3

Scan your own gene - just drag & drop

Save your sequence as a .fasta file, then DRAG it onto the SCAN-a-gene icon.
(Mac: double-click SCAN-a-gene.command, then pick your file.)



4

Read the result - worst first

You get a plain list. Fix the SERIOUS items first; each shows where it is and a suggested fix. Re-scan until they are gone.

SERIOUS - fix first
* Hidden lab cut-site
* Premature stop signal
Fix: change GGT->GGC ...

Needs Python 3.8+ (free): python.org/downloads - the launcher tells you if it is missing.

Scope: checks the mRNA 'message' only. Cannot see protein stability. No expression score. Chlamydomonas nucleus only.
A clean report is not a promise of expression. Full details in [USER_GUIDE.pdf](#).